

BinaxNOW[™] RSV Card

Fast, accurate detection to help reduce
the spread of RSV infections




BinaxNOW™ RSV Card



FEATURES AND BENEFITS

- **Results:** 15 minutes
- **Onboard extraction:** no extra steps
- **Sample types:** nasopharyngeal swab and nasal wash
- **Shelf life:** up to 24 months
- **One step:** no additional reagents needed

EASY TO USE



Refer to Package Insert for full instructions.

BENEFITS OF RAPID DIAGNOSIS

- Guide appropriate therapy
- Avoid unnecessary antibiotic use
- Reduce transmissions to others

WHY TEST FOR RSV?

- More than 50,000 hospital admissions for children per year nationally¹
- Clinically difficult to differentiate RSV from other serious respiratory viruses²

RSV COMPARISON GRID			ORDER INFORMATION	
	BinaxNOW™ RSV ³	QUIDEL® QUICKVUE® RSV ⁴	BinaxNOW™ RSV – CLIA WAIVED	
Format	Card	Dipstick	22 Test Kit	Order Number: 430-122
CLIA Waived	Yes	Yes	10 Test Kit	Order Number: 430-100
Sample Types	NP swab or nasal wash	NP swab, NP/nasal wash or NP aspirate	BinaxNOW™ RSV – MODERATELY COMPLEX	
Reagents	0	1	42 Test Kit	Order Number: 430-000
Timed Steps	1	2	NP SWAB ACCESSORY PACK	
Tests per Kit	22	20	22 swabs and 22 vials of elution solution	Order Number: 400-065
Kit Storage	Room temperature	Room temperature		
Shelf Life	24 months	24 months		
QC Requirements	Once per shipment	Every new operator, each new lot or every shipment (provided that each different lot in the shipment is tested)		

YOU TEST FOR ONE. WHY NOT BOTH?



TO LEARN MORE ABOUT RAPID DETECTION OF RSV, CONTACT YOUR LOCAL ACCOUNT REPRESENTATIVE, CALL 1-877-441-7440 OR VISIT ABBOTT.COM/POCT

NP = nasopharyngeal
RSV = respiratory syncytial virus

1. Midgley CM, Haynes AK, Baumgardner JL, et al. Determining the seasonality of respiratory syncytial virus in the United States: the impact of increased molecular testing. *The Journal of Infectious Diseases*. 2017;216(3):345-355.
2. Falsey AR, Walsh E. Respiratory syncytial virus infection in adults. *Clinical Microbiology Reviews*. 2000;13(3):371-384.
3. BinaxNOW™ RSV Package Insert. IN430002.
4. Quidel® QuickVue® Package Insert.

© 2022 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. COL-01161-03 08/22

