


Reservoirs



Nest reservoirs have the standard SBS size, which is suitable for all kinds of automatic liquid handling systems. The built-in groove at the bottom of the reservoirs can reduce the dead volume of reagent and ensure the consistency of the liquid volume being handled.


Features

- SBS International Standard
- Medical grade polypropylene (PP).
- Autoclavable at 121 °C/ 15psi.
- Non - pyrogen, DNase/RNase Free.
- Lot traceability.
- Low and Medium profile designs to satisfy your liquid handling needs.
- They come without caps, 10/pk, 50/cs.

Multi Well, 12 Channel, Low Profile

Description		Channel	Well Volume	Reservoir Height	Non-sterile Cat. No.	Sterile Cat. No.
Multi Well	Low Profile 31.4 mm Height	8	22 mL	31.4 mm	360101	360161
		12	15 mL	31.4 mm	360102	360162
	Medium Profile 44.5 mm Height	8	32 mL	44.5 mm	360201	360261
		12	22 mL	44.5 mm	360202	360262
Single Well	Low Profile 31.4 mm Height	96	195 mL	31.4 mm	360103	360163
		384	195 mL	31.4 mm	360104	360164
	Medium Profile 44.5 mm Height	96	290 mL	44.5 mm	360203	360263
		384	290 mL	44.5 mm	360204	360264
		8	290 mL	44.5 mm	360205	360265
		12	290 mL	44.5 mm	360206	360266


Reservoirs



127.76*85.48*31.4 mm

Multi Well, 8 Channel, Low Profile


Non-sterile Cat. No.	360101
Sterile Cat. No.	360161
Well Volume	22 mL
Packing	5/pk, 50/cs



127.76*85.48*31.4 mm

Multi Well, 12 Channel, Low Profile


Non-sterile Cat. No.	360102
Sterile Cat. No.	360162
Well Volume	15 mL
Packing	5/pk, 50/cs



127.76*85.48*44.5 mm

Multi Well, 8 Channel, Medium Profile


Non-sterile Cat. No.	360201
Sterile Cat. No.	360261
Well Volume	32 mL
Packing	5/pk, 50/cs



127.76*85.48*44.5 mm

Multi Well, 12 Channel, Medium Profile

Non-sterile Cat. No.	360202
Sterile Cat. No.	360262
Well Volume	22 mL
Packing	5/pk, 50/cs




127.76*85.48*31.4 mm

Single Well, 96 Troughs, Low Profile

Non-sterile Cat. No.	360103
Sterile Cat. No.	360163
Well Volume	195 mL
Packing	5/pk, 50/cs


Reservoirs



127.76*85.48*31.4 mm

Multi Well, 384 Troughs, Low Profile


Non-sterile Cat. No.	360104
Sterile Cat. No.	360164
Well Volume	195 mL
Packing	5/pk, 50/cs



127.76*85.48*44.5 mm

Multi Well, 96 Troughs, Medium Profile


Non-sterile Cat. No.	360203
Sterile Cat. No.	360263
Well Volume	290 mL
Packing	5/pk, 50/cs



127.76*85.48*44.5 mm

Multi Well, 384 Troughs, Medium Profile


Non-sterile Cat. No.	360204
Sterile Cat. No.	360264
Well Volume	290 mL
Packing	5/pk, 50/cs



127.76*85.48*44.5 mm

Multi Well, 8 Troughs, Medium Profile

Non-sterile Cat. No.	360205
Sterile Cat. No.	360265
Well Volume	290 mL
Packing	5/pk, 50/cs



127.76*85.48*44.5 mm

Single Well, 12 Troughs, Medium Profile

Non-sterile Cat. No.	360206
Sterile Cat. No.	360266
Well Volume	290 mL
Packing	5/pk, 50/cs