



BACTRACK® ELEMENT

BREATHALYZER


ESTIMATES BLOOD ALCOHOL CONTENT
IN SECONDS

ENGINEERED FOR MAXIMUM ACCURACY
AND PORTABILITY – IT'S BACTRACK'S
SMALLEST FUEL CELL BREATHALYZER EVER!

Perfect for

- Testing yourself, friends, and family
- Businesses, schools, hospitals, and law enforcement
- Learning how drinking alcohol affects your BAC
- Staying safe and consuming alcohol responsibly

ADVANCED FEATURES Automatic breath flow check and temperature check help ensure proper use.



TRUSTED BY PROFESSIONALS Uses police-grade sensor technology for trusted accuracy and consistency.

FAST RESULTS Blow into the mouthpiece and results are displayed in under 10 seconds.

PORTABLE DESIGN Under five inches tall and less than an inch thick – take it anywhere alcohol is consumed.



FOR MORE INFORMATION

BACTRACK® ELEMENT

FUEL CELL TECHNOLOGY IN A PORTABLE DESIGN




The BACTRACK Element – Sleek, Portable Design

The new BACTRACK Element quickly and easily tests for the presence of alcohol and provides a Blood Alcohol Content (BAC) estimate. Simply press power, blow through one of the included mouthpieces, and get results in seconds. The Element uses professional Xtend™ fuel cell sensor technology and an internal breath pump to provide trusted results. The Element's small size and affordability make it perfect for anyone – from individuals and businesses to professional groups and even law enforcement.

PRODUCT SPECIFICATIONS

| | |
|-----------------------------|--|
| Dimensions..... | 4.75 x 1.90 x 0.75 inches (12.1 x 4.8 x 1.9 cm) |
| Weight | 3.3 oz (94 g) With mouthpiece and batteries |
| Sensor Technology..... | Xtend™ fuel cell sensor technology |
| Detection Range | 0.00 – 0.40 %BAC |
| Mouthpieces | 2 included |
| Test Count | Displays total number of tests performed |
| Power Supply | Two AAA alkaline batteries, included |
| Battery Life | Approximately 1000 tests |
| Warm Up Time | 10-20 seconds, depending on last use |
| Response Time..... | 10 seconds |
| Sensor Accuracy | +/- 0.005 %BAC @ 0.1 %BAC |
| Operating Temperature | 32-104 °F (0-40 °C) |



- A** Breath Intake Area – versatile design allows quick use with replaceable, hygienic mouthpieces.
- B** Crisp 4-digit multifunction display – Easy to read test results and alerts.
- C** One Button Operation – Simply turn it on and begin testing in seconds.
- D** Compact size and sleek design - Easily fits into a pocket or purse for on the go lifestyles.

FOR MORE INFORMATION