


90x53 mm Urine Containers

- Sealed Lid, Stackable



Patient Label

Temp.Strip Attached


Label sealed
95Kpa compliant
Wide mouth size fit dip cards
Sealed, sterilized, Individually wrapped

Cat#	Descriptions	Package
JYB-SPUC-90mL	90mL Cups, with temp.strip & ID label, Sterilized	1pc/pk, 400/cs

100-120mL Vacuum Containers

Frosted labeling area
Integrated transfer device
Available in Temp.Strip for option
With needle built-in, OTC use
Easy to store and transport, safe and reliable.

- Vacutainers



Plain Tubes



Cat#	Descriptions	Package
JYB-EGCT1026	120mL vacutainers, yellow cap & sterile	1pc/pk, 250/cs
JYB-VAC-CUP-T	120mL vacutainers, yellow cap / temp.strip , sterile	1pc/pk, 250/cs
JYB-ST-U10P-100	10ml vacuum tubes, 16*100mm / non-additive, sterile	100pc/pk, 1200/cs
JYB-ST-U5P-100	5ml vacuum tubes, 13*100mm / non-additive, sterile	100pc/pk, 1800/cs
JYB-ST-U5P-75	5ml vacuum tubes, 13*75mm / non-additive, sterile	100pc/pk, 1800/cs

Distributed By:



2721 Loker Ave. W, Carlsbad, CA92010



<https://cliawaived.com/>



1-888-882-7739



info@cliawaived.com




Catalog No.	JYB-SPUC-90mL
Description	90mL Cups, with temp.strip & ID label, Sterilized
Built-in Needles (Cup)	No
Components	3oz/90ml urine cups / temperature strip / ID Patient Label
Package	Cups: individual wrapped, 400pcs/cs, 56*42*52cm, 7/8kgs/case
Frosted Area	No
Material	Polypropylene (PP)
95 Kpa Compliant	Yes
Temperature Strip	Yes
Sterile	EO/ Gamma
Label Custom	Yes
MOQ	400pcs
Color	White cap, clear cup / custom available
Certificate	COA/FDA/CE
validity Period	5 years
Sample available	Yes




Catalog No.	JYB-EGCT1026
Description	120ml vacuum urine cup, w/needle, graduation
Built-in Needles (Cup)	Yes
Components	4oz vacuum cups
Package	individual wrapped, 250pcs/cs, 64*31*46cm, 7.3kgs/case
Frosted Area	Yes, Name/Room No/Doctor/Date/Time.
Material	Polypropylene (PP)
Vacuum Pressure	Setting based on buyer's local area
Temperature Strip	No
Sterile	EO/ Gamma
Label Custom	Yes
MOQ	250pcs
Color	Yellow / custom available
Certificate	COA/FDA/CE
validity Period	5 years
Sample available	Yes




Catalog No.	JYB-VAC-CUP-T
Description	120ml vacuum urine cup, w/needle, graduation
Built-in Needles (Cup)	Yes
Components	4oz vacuum cups
Package	individual wrapped, 250pcs/cs, 64*31*46cm, 7.3kgs/case
Frosted Area	Yes, Name/Room No/Doctor/Date/Time.
Material	Polypropylene (PP)
Vacuum Pressure	Setting based on buyer's local area
Temperature Strip	Yes
Sterile	EO/ Gamma
Label Custom	Yes
MOQ	250pcs
Color	Yellow / custom available
Certificate	COA/FDA/CE
validity Period	5 years
Sample available	Yes



Catalog No.	JYB-ST-U10P-100
Description	10ml Vacuum Tube, 16*100mm, cap yellow, non-additive
Additive	Non, Plain Tube
Size	16*100mm
Package	100pcs/tray, 1200pcs/case. 51*41*26cm, 13/14kgs/case
Volume	10ml
Material	PET
Usage	For used with vacuum urine cup
Cap	Push Top type, Butly rubber inside
Sterile	EO/ Gamma
Label Custom	Yes
MOQ	100pcs
Color	Yellow / custom available
Certificate	COA/CE/ISO
validity Period	2 years
Sample available	Yes



Catalog No.	JYB-ST-U5P-100
Description	5ml Vacuum Tube, 13*100mm, cap yellow, non-additive
Additive	Non, Plain Tube
Size	13*100mm
Package	100pcs/tray, 1800pcs/case. 54*32*37cm, 14/15kgs/case
Volume	5ml
Material	PET
Usage	For used with vacuum urine cup
Cap	Push Top type, Butly rubber inside
Sterile	EO/ Gamma
Label Custom	Yes
MOQ	100pcs
Color	Yellow / custom available
Certificate	COA/CE/ISO
validity Period	2 years
Sample available	Yes



Catalog No.	JYB-ST-U5P-75
Description	5ml Vacuum Tube, 13*75mm, cap yellow, non-additive
Additive	Non, Plain Tube
Size	13*75mm
Package	100pcs/tray, 1800pcs/case. 45*36*29cm, 12/13kgs/case
Volume	5ml
Material	PET
Usage	For used with vacuum urine cup
Cap	Push Top type, Butly rubber inside
Sterile	EO/ Gamma
Label Custom	Yes
MOQ	100pcs
Color	Yellow / custom available
Certificate	COA/CE/ISO
validity Period	2 years
Sample available	Yes