

Nano-Check™ RSV Test


CLIA 
WAIVED

FDA
CLEARED

- ✓ FDA 510(k) Cleared (k240280)
- ✓ FDA CLIA Waived (CW240005)
- ✓ Anterior Nasal Swab
- ✓ Infant/Child: 6 months - 6 years old
- ✓ Adult Senior: ≥ 60 years old
- ✓ Results in 15 minutes
- ✓ No Extra Steps
- ✓ Relative Sensitivity: 83.8%
- ✓ Relative Specificity: 99.6%
- ✓ Made in USA



► Easy/Simple To Use



Contents

- 25 Test Devices
- 25 Reagent tubes with buffer
- 25 Sample collection swabs
- 1 Positive control swab
- 1 Negative control swab
- 1 Instructions for Use
- 1 Quick Reference Instruction

► Performance Characteristics

RSV subtype	Strain	LoD
A	Long	1.6×10^3 TCID 50/mL
	A-2	2.0×10^4 TCID 50/mL
B	18537	4.0×10^3 TCID 50/mL
	B1	3.5×10^3 TCID 50/mL

Refer to IFU for additional LOD / Analytical Reactivity information

Ordering Information

Cat. No.: ND-MD8152 (25 Test)

► Distributed By:

CLIAwaived.com

Tel: 1.858.481.5031

sales@cliawaived.com