



## STATEMENT

To whom it may concern,

We hereby declare about the minimum amount of Multi-Drug Rapid Test Cup. Based on the sample volume flex study, the product could yield correct and valid results when using a sample volume of 5mL or more for the cup application. Therefore, it is confirmed that the minimum amount of Multi-Drug Rapid Test Cup is 5mL.

We would also like to thank you for your support for making us the market leader in the sector of the IVD kit market and apologize for any shortcomings that may have occurred. We will always strive to do better. Stay safe.

Yours Sincerely,

A handwritten signature in black ink, appearing to read 'Darren Wang', written over a horizontal line.

Darren Wang  
Rapid diagnosis R&D Manager