

Are you frustrated with your Current Rapid Lateral Flow Fentanyl Test?

Did you know that Rapid Visual Read Lateral Flow Fentanyl Tests are only 50% Accurate and have a High level of Cross Reactivity?

Ultra Sensitive
1 ng/ml



BioSieve™ Toxismart FIA Reader with BioSieve™ Fentanyl FIA Pro Tests provide healthcare professionals with a fast and accurate way to detect fentanyl in urine samples. It is a portable fluorescence instrument for *in vitro* diagnostic use only.

Sample Concentration	Test Result	
	BioSieve™ Toxismart	Brand "A"
0.5 ng/mL	10(-) / 0(+)	10(-) / 0(+)
1 ng/mL	0(-) / 10(+)	10(-) / 0(+)
1.5 ng/mL	0(-) / 10(+)	10(-) / 0(+)
2 ng/mL	0(-) / 10(+)	7(-) / 3(+)
3 ng/mL	0(-) / 10(+)	4(-) / 6(+)
4 ng/mL	0(-) / 10(+)	2(-) / 8(+)
5 ng/mL	0(-) / 10(+)	0(-) / 10(+)

* Data on File

Test Result (n=365)	BioSieve™ Toxismart
Positive	0
Negative	365
*CLINICAL Specificity	100%

Studies show that BioSieve ToxiSmart FIA Pro System has a True Sensitivity of **1 ng/ml**.

Competitor Brand "A" has a Sensitivity between **3-5 ng/ml**.

BioSieve™ Toxismart FIA PRO System



Key Features

- Ultra-Sensitive 1 ng/ml
- Overall Accuracy: >99%
- Small Sample (3 Drops) of urine
- Portable: Use in Multiple Exam Rooms
- Easy to Use-minimal training required
- Handheld: light weight (5.76oz)
- Provides Printable Results
- Reimbursable
- Suggested Reimbursement CPT Code 80306-QW

Ordering information

Item code (REF)	Description	Size
VCUS-FIARDR	BioSieve™ Toxismart™ FIA Reader and BioSieve™ Toxismart™ Bluetooth Printer with charging cables, power adapter, SQC test device, user's manual and QRGs	1 Reader/Box
VCUS-FIARDRPROMO	BioSieve™ Toxismart™ FIA Reader and BioSieve™ Toxismart™ Bluetooth Printer with manuals, SQC test device, charging cables and 100 x BioSieve™ Fentanyl FIA Pro Tests	1 Promo/Box
VC-FYL1	BioSieve™ Fentanyl FIA Pro Test Kit, each box contains 25 test devices, 25 disposable pipettes, package insert and QRG	25 Tests/Kit
VC-FYCTL	BioSieve™ Fentanyl Controls, 5 ml Negative and 5 ml Positive Controls (1 Set)	1 Set/Box
VC-TPPR10	BioSieve™ Thermal Paper for use with BioSieve™ Toxismart™ Bluetooth Printer	10 Rolls/Box



VivaChek Biotech (Hangzhou) Co., Ltd.
Level 2, Block 2, 146 East Chaofeng Rd., Yuhang Economy
Development Zone, Hangzhou, 311100, Zhejiang, P.R. China
VivaChek Diagnostics, Inc.
1401 21st St STE R Sacramento, California, 95811
United States of America
TEL: (888) 860-3894
E-mail: info@vivachek.com
Website: <https://vivachekdx.com>

Number: 1606017206
Effective Date: 2024-11-22
©2024 VivaChek Biotech (Hangzhou) Co., Ltd.