



Advin Biotech, Inc.

10140 Mesa Rim Rd, Unit 1

San Diego, CA 92121, USA

Phone/Fax 1-858-866-8382 • www.advinbio.com •

info@advinbio.com

Mar 24, 2025

Re: BUP detects both the parent compound and metabolite

This letter clarifies the different cross-reactivity listings for BUP in package inserts 36001-CWII and G36002-CW.

The two package inserts, although both accurate, reflect different sets of cross-reactivity data based on the timeline of their creation and intended use.

The first insert, created in 2013 and submitted to the FDA, includes a single cross-reactivity. Due to its submission and approval process, making changes to this data is not straightforward.

The second insert, developed later for forensic use items, lists four cross-reactivities.

Despite these differences, both inserts use the same formula and chemicals, ensuring consistency in the product's composition.

Thank you for your business.

Quality Assurance

Advin Biotech