

The FastPack® IP hCG Immunoassay



FAST / SIMPLE / ACCURATE

The FastPack® IP hCG Immunoassay

The FastPack® IP hCG Immunoassay, for use on the FastPack® IP System, allows any testing facility to produce point-of-care results in 12 minutes with less than one minute of hands-on time.

- Enhanced quality of service for patients
- Ability to monitor and treat in one visit
- Improved patient compliance and outcome
- Laboratory accuracy and precision



SEKISUI
DIAGNOSTICS

Because every result matters™

The FastPack® IP hCG Immunoassay

Simple Process



1
*Dispense sample into the
FastPack® Injection Port*



2
*Place FastPack® into the
FastPack® IP Instrument*



3
Press the START button

The FastPack® IP hCG Assay Specifications:

- Quantitative measurement of hCG
- Lower limit: 1.8 mIU/mL
- Upper limit: 1000 mIU/mL; up to 10,000,000 with dilution
- Assay time: 12 minutes
- CPT CODE: 84702

ORDERING INFORMATION	CONFIGURATION	CATALOG NUMBER
FastPack® IP hCG	30 Tests	25000044
FastPack® IP hCG Calibration	Kit	25000035
FastPack® IP hCG Verification	Kit	25000033
FastPack® IP hCG Diluent	Kit/25	25000036
FastPack® Control (US)	Kit	25000056
FastPack® Control (International)	Kit	25000065

THE AMERICAS

Sekisui Diagnostics, LLC
4 Hartwell Place
Lexington, MA 02421
Phone: 800 332 1042
Fax: 800 762 6311
Email: questions@sekisui-dx.com

INTERNATIONAL

Sekisui Diagnostics (UK) Limited
Liphook Way, Allington
Maidstone, Kent ME16 0LQ, UK
Phone: +44 1622 607800
Fax: +44 1622 607801
Email: info@sekisui-dx.com

www.sekisuidiagnostics.com

SEKISUI
DIAGNOSTICS
Because every result matters™