


Upgrade your point-of-care hemoglobin testing with the HemoCue® Hb 801 System

	 HemoCue® Hb 201+ System	 HemoCue® Hb 801 System
Speed	Results can take up to 60 seconds	Results in less than 1 second
Microcuvettes	Reagent microcuvettes	Reagent-free microcuvettes for improved open vial stability
Microcuvette shelf life	18 months	30 months
Microcuvette stability	3 months for open vials	30 months for open vials *From date of manufacturing
Display	Small display, error codes for maintenance	Large display with intuitive on-screen symbols
Connectivity	None	Bluetooth and middleware option *Not approved for hospital use

Accuracy always matters

Hemoglobin testing is often done as a part of a routine doctor's visit or when a patient presents symptoms of anemia.¹ With a small drop of blood (~10 µL), the HemoCue® Hb 801 System determines an individual's hemoglobin value in less than one second. Calibrated against the International Council for Standardization in Hematology (ICSH) reference method, the HemoCue® Hb 801 System provides lab-quality results in CLIA-waived settings.



Each point on the diagram represents an actual hemoglobin test result. The closer to the line, the more closely related the result is to the "gold standard" reference method.²




Committed to accuracy since day one



For more than four decades we've been advancing what's possible at the point of care. Today, HemoCue is one of the most widely used point-of-care Hb systems in the world.³

Our caring commitment to our physicians, partners and patients remains as strong as it was on the first day of our journey.



- **One-Touch Resolution**
Get answers from our live customer service and technical support teams located at HemoCue America in California
- **In-Field or Virtual Service**
Get extensive support from your HemoCue representative virtually or right in your office
- **Online Education**
Easily stay up to date with comprehensive eLearning modules, competency reviews, and more

Product Order Number

121920

HemoCue® Hb 801 Analyzer

Make the HemoCue® Hb 801 System part of your workflow

Contact your HemoCue or preferred distribution representative

Learn more at hemocue.com/us/products/hematology/hemocue-hb-801-system/

1. <https://www.mayoclinic.org/tests-procedures/hemoglobin-test/about/pac-20385075>
2. HemoCue® Hb 801 System V&V Graphs, Figure 2; HCAB:000020941[2]
3. Sanchis-Gomar F, Cortell-Ballester J, Pareja-Galeano H, Banfi G, Lippi G. Hemoglobin point-of-care testing: the HemoCue system. J Lab Autom. 2013 Jun;18(3):198-205. doi: 10.1177/2211068212457560. Epub 2012 Sep 6. PMID: 22961038