

CDS 3PD Hematology Control

KIT LOT: 0812-502 EXP: 27-Apr-09
 Open Vial Stability: 14 days

LOW CONTROL		MINDRAY	
Lot 0812-011		3000 Series	
Parameters		Mean	± Range
WBC	$\times 10^9/L$	3.5	0.3
Lymph#	$\times 10^9/L$	1.1	0.5
Mid#	$\times 10^9/L$	0.3	0.3
Gran#	$\times 10^9/L$	2.2	0.8
Lymph%	%	31.1	5.0
Mid%	%	7.3	4.0
Gran%	%	61.6	5.0
HGB	g/dL	6.4	0.3
RBC	$\times 10^{12}/L$	2.53	0.12
HCT	%	18.5	2.0
MCV	fL	73.2	4.0
MCH	pg	25.3	1.5
MCHC	g/dL	34.6	2.2
RDW	%	14.8	3.0
PLT	$\times 10^9/L$	94	15
MPV	fL	10.1	1.5
NORMAL CONTROL		MINDRAY	
Lot 0812-012		3000 Series	
Parameters			± Range
WBC	$\times 10^9/L$	8.7	0.6
Lymph#	$\times 10^9/L$	3.7	0.6
Mid#	$\times 10^9/L$	0.5	0.5
Gran#	$\times 10^9/L$	4.4	1.0
Lymph%	%	42.7	5.0
Mid%	%	6.2	4.0
Gran%	%	51.1	5.0
HGB	g/dL	12.6	0.4
RBC	$\times 10^{12}/L$	4.23	0.18
HCT	%	36.5	2.5
MCV	fL	86.2	4.0
MCH	pg	29.8	1.5
MCHC	g/dL	34.6	2.2
RDW	%	13.3	3.0
PLT	$\times 10^9/L$	254	30
MPV	fL	9.9	1.5
HIGH CONTROL		MINDRAY	
Lot 0812-013		3000 Series	
Parameters			± Range
WBC	$\times 10^9/L$	20.2	1.8
Lymph#	$\times 10^9/L$	10.6	1.5
Mid#	$\times 10^9/L$	1.8	1.0
Gran#	$\times 10^9/L$	7.8	2.5
Lymph%	%	52.7	5.0
Mid%	%	8.8	4.0
Gran%	%	38.5	5.0
HGB	g/dL	17.7	0.5
RBC	$\times 10^{12}/L$	5.31	0.22
HCT	%	51.7	3.0
MCV	fL	97.4	4.0
MCH	pg	33.3	1.5
MCHC	g/dL	34.2	2.2
RDW	%	13.0	3.0
PLT	$\times 10^9/L$	583	60
MPV	fL	10.1	1.5



Distributed by: **CLIAwaived.com™**
 San Diego, CA 92121
 tel 858-481-5031
 toll free 888-882-7739
 www.cliawaived.com

