

RAPID DIAGNOSTICS

osom[®]

Trichomonas Rapid Test

FOR THE DETECTION OF *TRICHOMONAS VAGINALIS* ANTIGENS



www.osomtests.com



FEATURES:

- ✓ **ACCURATE;** 83% sensitivity vs. culture; 95% agreement with culture and wet mount combined
- ✓ **FAST;** Results in 10 minutes or less
- ✓ **SAMPLE;** Vaginal Swabs
- ✓ **OBJECTIVE;** Detects the antigen; removes uncertainty associated with live organism instability
- ✓ **TRAINING;** Online training and certification; www.osomtraining.com
- ✓ **CUSTOM CARE;** Technical support staffed with certified medical technologists

RAPID DIAGNOSTICS

www.osomtests.com


osom[®]

Trichomonas Rapid Test

Sekisui Diagnostics, LLC (AMERICAS)
4 Hartwell Place Lexington, MA 02421
Phone: 800 332 1042
Fax: 800 762 6311
Email: osomrapidtests@sekisui-dx.com
ISO 13485 Certified

Sekisui Diagnostics (UK) Limited (INTERNATIONAL)
50 Gibson Drive, Kings Hill, West Malling, Kent, ME19 4AF, UK
Phone: +44 (0) 1732 220022 Fax: +44 (0) 1732 220024
Email: info@sekisui-dx.com
ISO 9001 Certified, ISO 13485 Certified

TEST PROCEDURE / RESULTS



Insert swab and mix.



Let stand 1 MINUTE.



Insert test stick into extracted sample.



Read results at 10 MINUTES.



DESCRIPTION

OSOM[®] TRICHOMONAS RAPID TEST
TRICHOMONAS POSITIVE CONTROL KIT

OSOM[®] VAGINITIS TRANSPORT SYSTEM (TRICH & BV)

CONFIGURATION

25 TESTS

50 SWAB SETS

CATALOG

181 OR 181E (CE MARKED)
182 (CE MARKED)

185

STORAGE

ROOM TEMP
ROOM TEMP

ROOM TEMP

CPT CODE

87808QW

