

# HbA1c results in minutes

Informed conversations at the point of care



## A1CNow+®

The A1CNow+® system from Chek Diagnostics provides healthcare professionals with a fast and easy way of getting HbA1c results—allowing informed conversations about how patients are managing their diabetes in minutes, not days.

For more information, visit [www.a1cnow.com](http://www.a1cnow.com).

### Key Features

- National Glycohemoglobin Standardization Program certified
- CLIA waived
- Reimbursable

DISTRIBUTED BY:

**CLIA** *waived*™ Inc.

To place your order, call 888-882-7739.

# HbA1c results in minutes

Informed conversations at the point of care



## A1CNow+®

### More efficient than the lab

- Results in 5 minutes
- Small (5µL) blood sample
- Portable – use in multiple exam rooms
- Easy to use – minimal training required
- Reimbursable



A1CNow+® system  
(10-count test kit)  
CHEK- 3024



A1CNow+® system  
(20-count test kit)  
CHEK-3021

To place your order, call 888.882.7739.

DISTRIBUTED BY:

**CLIAwaived™** Inc.

11578 Sorrento Valley Road #26, San Diego, CA 92121  
858.481.5031 | FAX 801.720.7568 | [cliawaived.com](http://cliawaived.com)