



Clearly different.
Clearly better.



Clearview® Strep A Exact II Dipstick

Clearview® Strep A Exact II Dipstick

Rapid technology
provides direct
detection of Group A
Streptococcal Antigen.



Clearly different

The Clearview® dipstick is a patented, innovative technology that provides a highly sensitive and specific test for Group A Strep.

- **Rapid Response:** Results in as little as 5 minutes.
- **Simplicity:** CLIA waived, easy-to-use dipstick test available.
- **Reliability:** Built-in controls ensure accuracy.
- **Flexibility:** Results while patient waits.

Clearly better

Clearview® Strep A Exact II provides physicians with fast, accurate results for an early diagnosis than can result in:

- Shorter course of illness.
- Reduced severity of symptoms.
- Reduced spread of bacteria to others.
- Reduced likelihood of complications.




Clearview® Strep A Exact II
30 Dipsticks
Product No. 92002
CLIA Waived

Clearview® Strep A Exact II Dipstick Comparison Grid

Product Name	Clearview® Strep A Exact II Dipstick	Quidel QuickVue® Dipstick Strep A	Genzyme OSOM® Ultra Strep A
Format	30 Dipsticks	50 Dipsticks	25 or 50 Dipsticks
CLIA Status	Waived	Waived	Waived
Sample	Swab	Swab	Swab
Sensitivity	97%	92%	92.6%
Specificity	95%	98%	92.8%
Procedural Steps	5	5	10
Time to Result	5 minutes	5 minutes	7 minutes
Shelf Life	24 months	18 months	12 months
Storage	Room Temperature	Room Temperature	Room Temperature

Clearview® Strep A Exact II Dipstick Procedure



www.alere.com

Contact your local Account Executive 877-441-7440 | alere.com