



Know Drug Test Cup™

When you need to know



DETECTS 15 DRUGS

AMP 500	mAMP500
BAR 200	MDMA 500
BUP 10	MTD 300
BZO 200	OPI 2000
COC 150	OXY 100
EtG 500	TRAM 200
FEN 50	THC 50

K2

JWH and AM-2201

+ ADULTERANTS (CRE, pH, OX)

- For workplace testing (Not FDA Cleared)
- Fast, dark test lines
- Non-negative results in five minutes
- Donor-friendly cup opening
- Results can be photocopied
- Optional tests for adulterants
- Remarkable pricing and quality

The Drug Screen Cup is intended for the detection of drugs of abuse and/or metabolites in human urine for workplace and insurance screening purposes only, excluding tests intended for Federal drug testing programs (SAMHSA, DOT, US Military).

The test provides a preliminary result only; presumptive positive results should be confirmed using an alternate chemical methodology (such as GC/MS, LC/MS, GC/MS/MS and LC/MS-MS) if donor doesn't acknowledge drug use or if your policies require.

Call your sales representative now at 1-888-882-7739