

Distributed by

CLIAwaived™ Inc.

Product Sheet

SURVILANCE SAFETY LANCETS

	SAFETY LANCET	COLOR	DEPTH	GAUGE	BLOOD FLOW
	SLN100	Yellow	1.0mm	G26	Low Flow 5-10 ul
	SLN200	Grey	1.8mm	G23	Low Flow 10-20 ul
	SLN240	Orange	2.2mm	G22	Med Flow 20-40 ul
	SLN300	Pink	2.8mm	G21	Med-High Flow 40-60 ul
	SLB200	Green	1.8mm	G19	High Flow 75-100 ul
	SLB250	Blue	2.3mm	G18	High Flow 150-200 ul