

STAT-Site® M^{Hgb} Products

Catalog No. 900900

STAT-Site® M^{Hgb} Meter

Includes STAT-Site® M^{Hgb} meter, operating manual, quick reference guide and battery (extra battery also included)



Catalog No. 901025

STAT-Site® M^{Hgb} Test Cards

Includes 100 test cards (4 packages of 25 each) and a lot-specific CODE key



Catalog No. 900905

STAT-Site® M^{Hgb} Battery Pack

Includes 5 batteries (approximately 5,000 tests)

Catalog No. 503000

STAT-Site® Hemoglobin Controls

6 x 1.5 mL (3 low and 3 high level)

STAT-Site® M^{Hgb} Specs

Methodology _____	Reflectance photometry
Principle _____	Azidemethemoglobin
Calibration _____	Automatic (CODE Key)
Sample Type _____	Whole blood
Sample Size _____	1 drop (approximately 12 µL)
Linearity _____	6 to 21 g / dL
Within Run Precision CV _____	3.5% to 4.9%
Total Precision CV _____	2.9% to 4.2%

Easy as one, two, three.

Step 1 INSERT THE TEST CARD

Step 2 APPLY THE SAMPLE

Step 3 READ THE RESULT



CPT CODE: 85018QW