




User Guide




Remove protective cap from the safety lancet.



Place red raised platform end of safety lancet on the test site.



Gently push the safety lancet down against test site to activate the lancet mechanism.



Once used, simply discard lancet into refuse container. Do not use lancet if damaged.

Why the SurgiLance™ Safety Lancet is Superior

Simple. No arming is required. Lancets are color-coded for easy recognition.

Safe. The needle or blade lancet is safely concealed before and after use.

Innovative. SurgiLance™ safety lancets makes blood sampling the easiest and safest method. Once the lancet is used, it is rendered inoperable, providing added safety for patient and clinician.

Efficient. The high-speed delivery and penetration method minimizes patient pain and clinician error.

Latex-free. SurgiLance™ safety lancets are hypo-allergenic.

www.medipurpose.com/surgilance

Americas

3883 Rogers Bridge Road NW
Suite 501
Duluth, GA 30097
Tel: +1 770 448 9493

Asia

15 Hoe Chiang Road #12-02
Tower Fifteen
Singapore 089316
Tel: +65 63451588

Europe

England, UK
Tel: +44 330 122 2630

MediPurpose™